

Adjustable Therapy Bench

CURRENT V1 BUILD GUIDE



IMPORTANT SAFETY NOTE

Confirm dimensions, positioning, hardware, supervision, and suitability with the child's therapist or care team. Inspect before every use. This guide is not medical advice or a substitute for commercial equipment when prescribed.



Materials, tools and cut list

TOOLS

- Circular saw or track saw
- Jigsaw
- Drill/driver
- Router or sanding tools
- Clamps and square
- Tape measure and pencil
- Eye and hearing protection

MATERIALS

- 3/4-inch birch plywood
- Bolts sized for leg channels
- Washers and threaded knobs
- Wood screws or pocket-hole screws
- Wood glue
- Sandpaper: 120 and 220 grit
- Optional non-slip top surface

CORE CUT LIST

PART	QTY	DIMENSIONS
Bench top	1	24 × 10 in
Inner legs	2	10 × 7 in
Outer legs	2	7 in high; 10 in top, 14 in bottom
Hardware	2 sets	Front and rear adjustment bolts

LAYOUT NOTES

- Keep the top grain running along the 24-inch length.
- Round all corners and exposed edges before assembly.
- Outer legs use a vertical channel with four horizontal locking slots.
- Front and rear bolt positions can differ to create tilt.
- Test the full adjustment range before finishing the wood.

Build sequence

1

Cut and shape the parts

Cut the top, two inner legs, and two trapezoid outer legs. Mark all pieces before drilling.

2

Create the adjustment channels

Lay out one vertical channel and four horizontal locking slots on each outer leg. Match both legs precisely.

3

Drill and smooth the slots

Drill starter holes, cut the channels with a jigsaw, then sand or route all slot edges smooth.

4

Attach the inner legs

Position the inner legs perpendicular to the top. Fasten securely using glue and suitable screws.

5

Fit the outer legs

Place each outer leg around its matching inner leg. Install bolts, washers, and knobs without binding.

6

Test every position

Check the lowest setting and all four slots. Confirm both sides lock fully and the bench does not rock.

7

Finish and inspect

Round edges, sand through 220 grit, and apply a child-appropriate finish if desired.



Final inspection and use

PRE-USE CHECKLIST

- ☐ All plywood edges are rounded and splinter-free.
- ☐ All bolts, washers, and knobs are installed correctly.
- ☐ Each locking slot fully captures the bolt.
- ☐ The bench sits level without rocking in every setting.
- ☐ No hardware protrudes into the user area.

ADAPT WITH GUIDANCE

A child's positioning, strength, range of motion, balance, and equipment needs are individual. Changes to height, angle, padding, restraints, or use should be reviewed with an appropriate professional. Stop using the bench immediately if any part loosens, cracks, shifts, or no longer fits the intended use.

SHARE WHAT YOU LEARN

Built For Them improves plans through real build feedback. Share photos, measurements, corrections, and lessons at hello@builtforthem.org.